

# Detection and quantification of rare mutations with massively parallel sequencing

Isaac Kinde, Jian Wu, Nick Papadopoulos, Kenneth W. Kinzler and Bert Vogelstein

Johns Hopkins Kimmel Cancer Center  
June 7, 2011


Tim Farrell  
BF831: Translational Bioinfo Seminar  
October 14, 2015

# Motivation

- Low-prevalence mutations are difficult to identify with confidence
  - In clinical setting, prevalence can be  $<0.01\%$
- Illumina, and sequencing technologies more generally, still have error rates high enough to mask these rare mutations
  - estimates ranged from 1% to 0.05% in 2011
  - $\sim 0.01\%$  in 2015 (1)
- Many aspects of sequencing factor into this error:
  - Library prep
  - Choice of primers
  - Base-calling
  - Read length

# Safe-Sequencing System (“Safe-SeqS”) Methodology

1. Assign unique id (UID; aka barcodes or indices) to each DNA template
2. Amplify of each uniquely-tagged template:
  - Creates many daughter molecules = “UID family”
3. Sequence
4. If 95% of any UID family has identical mutation → “supermutant”





# UID Assignment: Endogenous

Randomly sheared DNA  
inherently contains unique  
seq at either end → call UID

2 approaches:

- Capture: many genes simultaneously
- Inverse PCR: single gene in depth



# Endogenous UID Experiments

## Capture:

- 15000 normal cells
- 2.6kb from 6 genes
- SafeSeqS: error lowered 70 fold

## Inverse PCR:

- 1750 normal cells
- 1 gene of interest
- 1057 unique molecules assessed
- SafeSeqS: conventional results likely all FP

	Capture	Inverse PCR
Conventional analysis		
High-quality base pairs	106,958,863	1,041,346,645
Mean high-quality base pairs read depth	38,620x	2,085,600x
Mutations identified	25,563	234,352
Mutations/bp	2.4E-04	2.3E-04
Safe-SeqS analysis		
High-quality base pairs	106,958,863	1,041,346,645
Mean high-quality base pairs read depth	38,620x	2,085,600x
UID families	69,505	1,057
Average no. of members/UID family	40	21,688
Median no. of members/UID family	19	4
Supermutants identified	8	0
Supermutants/bp	3.5E-06	0.0


# UID Assignment: Exogenous

2 cycles of UID assignment during PCR

12-14 random nucleotides b/t forward primer and seq serves as UID

Amplification: 25 PCR cycles

Gives a uniquely-tagged, dsDNA fragment for each of the original template strands



# Exogenous UID Experiments

## CTNNB1:

- 3 individuals
- SafeSeqS: reduced apparent error 24 fold

## Mitochondrial DNA:

- 7 individuals
- SafeSeqS: 15 fold reduction

# Discussion

How do we know conventional analyses represent artifacts rather than true mutations?

- consistency of results
- polymerase fidelity experiment results consistent with previous biological assays

Limitations:

- Exogenous UID strategy cannot be applied to multiple amplicons from sample containing a limited # of templates
- Still does not overcome error introduced in amplification



# Future directions

# References

- (1) Schrimmer et al. 2015. Insight into biases and sequencing errors for amplicon sequencing with the Illumina MiSeq platform. *Nuc Acid Res.* doi: 10.1093/nar/gku1341.